

AmTryke®

Summer
2010



New and Improved
Equipment for
All Ages!

1-888-AMTRYKE
www.ambucs.org

Amtryke therapeutic tricycles featuring the combination of hand/foot motion are available in 3 sizes: AM-9, AM-12 and AM-16. Many accessories allow these trikes to adapt to a wide variety of sizes and ability.



AM-9S or AM-9XS can be converted to a hands only tryke with the addition of the foot platform. This will provide for riders that do not need the foot rotation component. *Coming Soon!*



AM-9S with AM-9 Foot Platform

AM-9 with Snappy Seat
#50-HFC-0102S with snappy seat
#50-HFC-0102XS with snappy seat

The Snappy Seat can be added to the AM-9S or AM-9XS hand/foot AmTryke. The adjustable Snappy Seat provides added stability for those children needing additional trunk support.



AM-9S With Snappy Seat

“ We use the hand-foot combination to improve coordination, strength and range of motion on a variety of subjects with unlimited diagnosis. The combination of using all four extremities helps with any weakness in any area. Children with generalized weakness can self propel despite the fact they are so very weak.” ~Sue Haywood, PT

AM-9 Hand / Foot AmTryke® Therapeutic Tricycle

AM-9S #50-HFC-0102S

AM-9XS #50-HFC-0102XS

This is the smallest tryke that features the combination hand/foot drive. It is designed to build strength and coordination for a wide range of children with disabilities. The hand and foot cycles are intended for early intervention by stimulating reciprocal movement and increasing range of motion. This tryke is suitable for children with a height of 28-36 inches.



AM-9S

- 4 inch hand crank arms
- Riders Height : 28-36 inches
- Arm Length (top of arm to middle of hand) 14-20 inches
- Leg Length (inseam to bottom of shoe) 14-20 inches

AM-9XS

- 2 ½ inch hand crank arms
- Forward Seat Extender Plate
- Riders Height: 24-36 inches
- Arm Length (top of arm to middle of hand) 10-20 inches
- Leg Length (inseam to bottom of shoe) 12- 20 inches

Standard Features

- Adjustable step-through frame
- Continuous chain mechanism for hand/foot movement
- Parental steering control on front of steering column
- Two-section adjustable frame; quick-turn knob allows frame to be quickly and easily locked into place
- AM-9 foot plates with heel trap, instep and toe straps
- Adjustable bucket seat with safety lap T-strap Two-section push bar with round plastic handle Safety helmet included.
- Maximum weight is 150lbs
- Assembly required



AM-9 footplate detail

Accessories:

Snappy Seat, Forward Extender Plate, 2 ½ inch crank arms, Hand Strap, Wrist Wraps, Wrist Brace Holding Mitts, Vertical Hand Grips, Brachial Plexus Kit, Foot Cups, Pedal Blocks, AM-9 Pedal Leveler ,Knee Adductor Positioning Strap



AM-9XS (Features XSMALL components) #50-HFC-0102XS

Xsmall components added to the AM-9 allow smaller children to ride. 2 ½ inch cranks replace the standard 4 inch cranks allowing a smaller range of motion for arms. A forward seat extender plate is also included to allow the bucket seat to be mounted closer to the front of the frame for shorter inseam/legs.

*AM-9XS Features the XSmall components
2 ½ crank and forward seat*



AM-12 and AM-16 saddle seat and plastic seatback with 2 chest straps

AM-12 and AM-16 Hand/Foot AmTryke Therapeutic Tricycle*

AM-12 #50-HFC-0210

AM-16 #50-HFC-0411

The AM-12 and AM-16 AmTrykes feature the combination hand/foot drive. Designed to build strength and coordination for a wide range of children with disabilities. The hand /foot cycles are intended to stimulate reciprocal movement and increase range of motion. Both trykes feature four-way adjustable seating that allows the seating option to go up/down and forward/back to adjust for a wide variety of rider sizes. Safety steering lock pin for parental control. Foot plates feature a heel trap, with toe and instep strap to secure riders feet to the pedals. Loading lock allows safe and easy transfer of riders. The AM-12 and AM-16 Model AmTrykes can be adapted with a wide variety of seating options and accessories so that most riders can feel safe and secure. *patent pending.



Seating can be positioned forward or back along the frame using large knobs



Seating height can be quickly adjusted with a spring loaded knob



Hinged Front End with 3 Positions for adjustable hand crank/ foot crank to seat alignment

Standard Features

- Continuous chain mechanism for hand/foot movement
- Safety Steering lock pin for graduated steering: straight, 20 degrees left/right, free steering
- Four-way adjustable seating- up/down and forward/back
- Unique hinged front end with three positions for adjustable hand crank/foot crank to seat alignment
- Self-righting pedals with heel trap, instep and toe straps
- Loading brake for easy transfer of rider
- Includes safety chest and lap straps
- Low center of gravity
- Available in red only
- Safety helmet included
- Dismantles for transport via car trunk, etc.
- Assembly required



AM-16 with saddle seat and 1400 seat back

AM-12

- 12 inch front and rear wheels
- Riders Height 36 – 42 inches
- Arm length (top of arm to middle of hand) 14-20 inches
- Leg length (inseam to bottom of shoe) 15-24 inches
- Maximum weight 150 lbs

AM-16

- 16 inch front and 12 inch rear wheels
- Riders Height 41 inches and greater
- Arm length (top of arm to middle of hand) 16-24 inches
- Leg length (inseam to bottom of shoe) 20-28 inches
- Maximum weight 175 lbs.

Accessories:

Snappy Seat (AM-12 only), Bench Seat, Tractor Seat, 1400 Seatback, 1400 Lumbar Pad, 1400 Head Rest, Full Padded Back, Hand Strap, Wrist Wraps, Wrist Brace Holding Mitts, Vertical Hand Grips, Brachial Plexus Kit, Foot Cups, Pedal Leveler Pulley, Exercise Pedals, Expanding Pedals, AM Foot Platform, Knee Separator, Knee Adductor Positioning Strap, Push Bar, Safety Flag



AM-12 with Snappy Seat



Safety Pin Lock for 3 Steering Positions: Straight, 20 Degrees right/left, free steering



AM-16 with bucket seat

Oceana from Kettler

#50-FC-0250

The Oceana is appropriate for children up to 36 inches in height. With a fixed front wheel, the pedals will always rotate which in turn will teach your child the motion of pedaling as well as building strong leg muscles. Standard tryke comes with AM-9 pedal plates with foot straps, a safety strap and helmet.

Accessories:

Snappy Seat, Forward Extender Plate, Hand Strap, Wrist Wraps, Wrist Brace Holding Mitts, Foot Cups, Pedal Blocks, AM-9 Pedal Leveler, Knee Adductor Positioning Strap



Oceana

Standard Features

- Quick-Adjust frame with telescoping multiple frame adjustments
- Parental control steering lock system locks handlebars into place
- Safety-Trax wheels with dual layer safety tread wheel pattern
- AM-9 pedal plates with heel trap, toe and instep straps
- Push bar
- Safety helmet
- Assembly required

Snappy

#50-FC-0150

The Snappy is suitable for children approximately 34 to 40 inches in height. Fixed front wheel, the pedals will always rotate which in turn will teach your child the motion of pedaling as well as building strong leg muscles. Adjustable seating allows positioning forward/back for multiple leg lengths. Snappy seat features padded bench seat and contoured back/head support which can be positioned forward/back. Low step through frame and loading brake make transferring riders safe and easy.

Standard Features

- 12- inch pneumatic front and rear wheels
- Safety Steering lock pin for graduated steering: straight, 20 degrees left/right, free steering
- Dismantles for transport via car trunk, etc.
- Adjustable seating- forward/back
- Pedal leveler pulley
- Low center of gravity
- Available in red only
- Safety helmet included
- Assembly required
- Maximum weight 150 lbs.
- Riders Height: 34- 40 inches
- Loading brake for easy transfer of rider
- Spring front wheel to bring steering back to center
- Leg Length (Inseam to bottom of shoe) 15-20 inches
- AM-9 pedal plate with heel trap, instep and toe straps
- Arm Length (Top of arm to middle of hand) 10-18 inches



Snappy with Snappy Seat



Oceana with Snappy Seat

Accessories:

AM-9 Bucket Seat, Hand Strap, Wrist Wraps, Wrist Brace Holding Mitts, Foot Cups, Pedal Blocks, Knee Adductor Positioning Strap, H-Harness, Lateral Supports, Safety Flag

The AmTryke Foot Cycles were developed in response to requests from therapists for a traditional tricycle for riders with special needs. Nine models of foot cycles- Oceana, Snappy, ProSeries 1412, 1416, 1420 and the 2700 Series Single and 7 Speed will provide many choices to provide the best match for the size and ability of a rider. These trykes can be adapted with a wide variety of seating options and accessories so that most riders can feel safe and secure.

AmTryke ProSeries 1412

#50-FC-1412

This positive drive tryke is designed for children approximately 36 to 42 inches in height and is built on 12-inch front and rear pneumatic wheels. Four way adjustable seat positions seating option up/down and front/back to accommodate different size riders. Quick release cam levers make handlebar positioning a snap. Loading lock allows safe transfer of riders on and off. Safety steering lock pin for parental control. The ProSeries 1412 can be adapted with a wide variety of seating options and accessories so that most riders can feel safe and secure.

Standard Features

- Positive drive action means foot crank is constantly in motion for full therapeutic effect.
- Four seat post positions
- Weighted self righting pedals
- Pedal plates with heel trap, instep and toe straps
- Loading lock for easy transfer of rider
- Saddle seat with lap strap and plastic back with 2 chest straps
- Quick release cam levers for easy handlebar adjustments
- Safety steering lock pin for three graduated steering positions: straight, 20-degrees left/right, free steering
- Safety helmet included
- Steel wire basket
- Assembly required
- Arm Length (Top of arm to middle of hand) 14-22 inches
- Leg Length (Inseam to bottom of shoe) 16-24 inches
- Riders Height: 36-42 inches
- Maximum weight 125 lbs



ProSeries 1412 with saddle seat and plastic back with 2 chest straps

Accessories:

Bucket Seat, 1400 Seatback, 1400 Lumbar Pad, Full Padded Back, H-Harness, 1400 Head Rest, Hand Strap, Wrist Wraps, Wrist Brace Holding Mitts, Foot Cups, Pedal Blocks, Exercise Pedals, Expanding Pedals, Pedal Leveler Pulley, Knee Adductor Positioning Strap, Knee Separator, Push Bar, Pull Bar, Safety Flag

ProSeries 1412 with Saddle seat and 1400 Seatback



ProSeries 1412 with bucket seat





ProSeries 1416 and ProSeries 1420

ProSeries 1416 #50-FC-1416
ProSeries 1420 #50-FC-1420

The ProSeries 1416 and ProSeries 1420 positive drive trykes provide full support to riders with low tone, trunk control. They can be adapted to meet the need of a wide size of riders, as four seat post positions accommodate multiple leg lengths. The saddle seat and 1400 seatback can be positioned up/down and in/out. The adjustable lateral supports with chest strap add additional trunk support. Quick release cam levers allow

positioning the telescoping handlebars to any height or distance to the rider. Expanding pedals with 2 straps adjust from 7 inches to 10 inches and the positioning over the pedal axis is also adjustable. A loading brake makes transfer of rider safe and easy. The ProSeries 1416 and 1420 can be adapted with a wide variety of seating options and accessories so that most riders can feel safe and secure.



1420 with wheel chair



Standard Features

- Positive drive action means foot crank is constantly in motion for full therapeutic effect.
- Safety steering lock pin for three graduated steering positions: straight, 20 degrees left/right, free steering.
- Four seat post positions to accommodate various rider sizes
- 1400 seat back (padded back rest and two adjustable lateral supports)
- Adjustable length expanding pedals with heel trap, toe and instep straps
- Foot support levelers
- Adjustable telescoping handlebars
- Loading lock brake for easy transfer of rider
- Detachable front frame for easy storage and transport
- Steel wire basket and helmet included
- Assembly required
- Optional Feature-Hand Brake



Loading Lock

ProSeries 1416

- 16 inch pneumatic front and rear wheels
- Riders Height: 41-60 inches
- Arm Length (top of arm to middle of hand) 16-24 inches
- Leg Length (inseam to bottom of shoe) 18-26 inches
- Maximum Weight 175 lbs

ProSeries 1420

- 20 inch pneumatic front and rear wheels
- Riders Height: 55 inches and greater
- Arm Length (top of arm to middle of hand) 18-26 inches
- Leg Length (inseam to bottom of shoe) 20-28 inches
- Maximum Weight 250 lbs



Expanding Pedals

Accessories:

Bucket seat (ProSeries 1416), Bench Seat, Tractor Seat, Wheel Chair Seat (ProSeries 1420) 1400 Lumbar Pad, Full Padded Back, H-harness, 1400 Headrest, Wrist Wrap, Wrist Brace Holding Mitt, Hi Rise Handlebars, MX Freestyle Handlebars, Footcups, Pedal Blocks, Exercise Pedals, Knee Separator, Knee Adductor Positioning Strap, Safety Flag



ProSeries with bucket seat



ProSeries with tractor seat

2700 Series AmTryke® Therapeutic Tricycles

#50-FC-2701 20 Inch Wheel Single Speed

#50-FC-2707 20 Inch Wheel 7 Speed

#50-FC-2721 24 Inch Wheel Single Speed

#50-FC-2727 24 Inch Wheel 7 Speed

The 2700 Series was developed to satisfy the requests of therapists who were seeking a cycle for higher functioning riders needing the full support provided by a wide wheelbase and stability of a tricycle. These tricycles are intended for use with teens/ adults with disabilities or to provide mobility and exercise opportunities for adults/ seniors. The low step through frame makes transfer safe and easy. The Safety steering lock pin allows control for steering: straight, 20 degrees right/left, or free steering. Seating options and accessories can adapt these trykes for a wide range of riders and abilities.

2701 20 Inch Wheel

2721 24 Inch Wheel

Both Available Single and 7 Speed



Standard Features

- Step-through frame for easy transfer on to tryke
- Wide saddle seat
- Steering lock pin for three graduated steering positions: straight, 20-degrees left/right, free steering
- Side pull caliber brake with built in lock for easy transfer of rider
- Features both front caliber brake and rear band brake
- Front and rear reflectors
- Large wire basket
- Bell
- Safety helmet included
- Assembly required

2700 Series 2701

- Single speed free wheel hub
- 20 inch front and rear wheels

2700 Series 2707

- 7 speed Shimano
- 20 inch front and rear wheels

2700 Series 2721

- Single speed free wheel hub
- 24 inch front and rear wheels

2700 Series 2727

- 7 speed Shimano
- 24 inch front and rear wheels



2700 Series with Tractor Seat



2700 Series with 2700 Seatback



2700 Series with 1400 Handle Bar



2700 Series with 2700 Seat Back/Laterals



Steering Pin Detail

Accessories:
Bench Seat, Tractor Seat, 2700 Seat Back/Laterals, Wrist Wraps, Wrist Brace Holding Mitt, 1400 Handlebars, Exercise Pedals, Expanding Pedals, 2700 Pedal Leveler Pulley

AmTryke Sizing Chart

Motion Controlled	Hand/Foot	Hand/Foot	Hand/Foot	Hand/Foot	Foot	Hand	Hand							
AmTryke Model	AM-9XS	AM-9S	AM-12	AM-16	Oceania	Snappy	1412	1416	1420	2701	2721	HP1000	1424	
Rider Inseam	12" to 20"	14" to 20"	15" to 24"	20" to 28"	14" to 20"	15" to 20"	16" to 24"	18" to 26"	20" to 28"	27" to 32"	28" to 36"	29" to 36"	Min. 23"	
Rider Arm length	9" to 20"	14" to 20"	14" to 20"	16" to 24"	14" to 20"	10" to 18"	14" to 22"	14" to 22"	18" to 26"			23" to 29"	23" to 36"	
Rider Weight (lbs)	150	150	150	175	150	150	125	175	250	275	275	250	250	
Tryke Weight (lbs)	33lbs	33lbs	45lbs	55lbs	33lbs	45lbs	72lbs	74lbs	74lbs	72lbs	74lbs	74lbs	74lbs	
Wheel Size	10"	10"	12"	16"	10"	12"	12"	16"	20"	20"	24"	Front 16" Rear 20"	20"	
Frame Height (Without Seat)	11"	11"	11	12.5"	11"	11"	13"	16.5	18.5	25	28	17	14.5	
Total Length (with push bar add 24")	26"	26"	39.5	45"	26"	36"	43"	53	64"	66"	72"	72"	64"	
Tryke Width	19.5"	19.5"	25	31"	19.5"	24"	27"	28	30"	29.5"	29.5"	29.5"	29.5"	



Coming Soon!
Swirlygig Hand Cycle
Kids Love it!

Helmet Sizing	
Style	Size in Inches
XXS	18.5" to 19.5"
XS	20.5" to 22"
S/M	22" to 23 5/8"
Lg/XL	23 5/5" to 25 3/4"

Introducing the AmTryke Hand Cycles

1424 Community Cruiser AmTryke® Therapeutic Tricycle

#50-HC-1424

The 1424 Community Cruiser provides a good mobility option for larger riders and offers an adjustable easy slide wheel chair seat for easy transfer of rider. Riders who are comfortable with the normal wheelchair position find this seat familiar and easy to adjust. The hand crank tryke operates easily, and includes a 3-speed Shimano hub for ease in riding on gentle slopes and hilly terrain. A tight turning radius offers the rider easy maneuverability. Adjustable foot platforms provide secure foot positioning. This community cruiser features a large basket for shopping, a coaster brake for additional safe riding and a front hand brake.

Standard Features

- 3 Speed Shimano Hub Control
- Front Caliper Brake
- Coaster Brake
- Auto Center of Steering
- 24 inch pneumatic front wheel
- 20 inch pneumatic rear wheels
- Adjustable Foot Platforms
- Easy Slide Wheel Chair Seat Adjustment
- Basket
- Helmet
- Assembly required
- Maximum weight 275 lbs

Accessories:

Tractor Seat, 2701 saddle seat, 2721 saddle seat, 1400 Seatback set with saddle seat, Safety Flag



1424 Community Cruiser



HP1000 Recreation Cycle

#50-FC-1000

The HP-1000 hand cycle is a full featured recreation cycle. It boasts a fully adjustable recumbent seat that will fit almost any adult rider. The cycle features a 7- speed Shimano hub paired with rear disk and front hand brakes. The cycle has adjustable leg rests. It is designed for comfortable riding on paved trails and roads. The multiple speeds will allow the rider to keep pace with family members and friends in most situations and achieve higher speeds when appropriate. The turning radius is wider than the 1424 Community Cruiser, at 10 feet.



Standard Features

- Rear disc brake
- Front caliper brake
- 7 Speed Shimano hub control for advanced riders
- Conforming nylon mesh seat
- Adjustable leg supports
- Quick release seat for easy adjustment
- Spring centered steering
- Helmet
- Assembly required
- Maximum weight 275 lbs.



AM-9 Hand Cycle

The AM-9S and AM-9Xs can be converted to a hand cycle with the addition of a foot platform.



AM-12 and AM-16 Hand Cycles

The AM-12 and AM-16 can be converted to a hand cycle with the addition of a foot platform.

Accessories

Accessories allow adaptations for the trykes so that most riders can ride the AmTryke® Therapeutic Tricycle. Accessories for seating options, trunk support, hand, and pedal options provide for a wide variety of riders and abilities.

Seating

Tractor Saddle Seat

#30-01-0201

The tractor saddle seat provides a large contoured area for stable seating. Fits the AM-12, AM-16, ProSeries 1416 ,1420 and the 2700 Series 20 inch and 24 inch model trykes.

Tractor saddle seat for the AM-12, AM-16, ProSeries 1416 or 1420 will need the following: #30-01-0212 AM/1400 Series Support Bracket

Tractor saddle seat for the 2700 20" or 2700 24" A rear support bar is needed.
#30-01-0202-19 (19" Length support bar for riders from 5'1" to 5'5")
#30-01-0202-23 (23" Length support bar for riders from 5'4" to 5'8")
#30-01-0202-26 (26" Length support bar for riders from 5'7" and up)



Snappy Seat

#30-04-1500

Fits: AM-9S, AM-9XS, AM-12 and is standard on the Snappy

The snappy seat provides trunk support for the smallest riders. The 2 part seating system's base is a 8 inch wide and 10 inch deep bench seat . The 2nd part of the seat is the contoured padded back and head support. The head support can be raised up and down. The padded back/head support can be easily adjusted at any position along the 10 inch seat. An H-Harness and two sizes of laterals can also be added to this adjustable seat to provide more stability.



Bucket Seat

#30-02-0320

Fits: AM-12, AM-16, ProSeries 1412 and 1416

The bucket seat comes with a 50 inch black hook and loop lap strap and 40 inch chest strap. The bucket seat is a popular choice for riders who need more support, but not a full trunk support seat back and laterals. The molded plastic seat measures 13 inches wide, 8 inches seat depth and 11 inch seat back. The bucket seat replaces the saddle seat.



Wheel Chair Seat

#30-03-0001

Fits: ProSeries 1420 and is standard on the 1424

The wheel chair seat provides for an extra roomy comfortable seat. The padded seat measures 18 inches wide and 16 inches deep. Safety belt included.



Saddle Seats

Variety of saddle seat sizes will allow a comfortable choice in seating.

#30-01-0109 AM/1400 Saddle Seat 11 inches wide by 9 inches deep, textured

#30-01-0110 2701 Saddle Seat 10.5 inches wide by 10 inches deep

#30-01-0111 2721 Saddle Seat 13 inches wide by 12 inches deep

Bench Seat

#30-05-0627

Fits: AM-12, AM-16, ProSeries 1416, 1420, 2700 Series 20 and 24 inch.

The bench seat is used when more hip support is needed than the standard saddle seat would provide. The bench seat measures 13 1/2 x 9 1/2 inches and is made of easy to clean polypropylene.

The bench can be tilted upwards for better hip stability.

The bench seat replaces the saddle seat. Weight limit 250lbs.



“This bike is going to help my daughter get stronger. It also gives her something in common with other children her age.”

~The Jensen Family

Trunk Support

1400 Seat Back

#30-13-1400

Fits: AM-12, AM-16, ProSeries 1412 and is standard on the ProSeries 1416 and 1420.

This versatile seat back uses a double tube support system. The 8 inch wide by 10 inch tall padded back can be adjusted by large knobs up/down the double tube. Quick release cam lock levers allow the 2 removable padded lateral supports adjustments up/down, forward/back and swing in/out. Plus a strap secures lateral supports in front of rider. Features built-in push grip for caregiver or parent. The 1400 seat back can also add the full padded seat back, lumbar pad, and 1400 head rest. The 1400 seat back can be used with a saddle seat, bench seat, or tractor seat.



Full Padded Back

#30-12-0335 Fits: 1400 Seat Back

The full padded back provides additional trunk support for riders with low muscle tone. Replaces the padded back that is standard on the 1400 seat back. Comes with brackets for mounting H-Harness. Measures 10 inches wide by 15 inches tall.

Coming soon, Full padded back for AM-12 and AM-16.



H-Harness

#30-11-0336 11.5 inch

#30-11-0337 8 inch

Provides additional trunk support. Fits full padded back.



1400 Lumbar Pad

#30-12-0334 Fits: 1400 Seat back

This 9 inch wide by 4 inch tall pad can be mounted on the 1400 seat back as additional support for the standard seat back pad. It can be positioned below the standard pad as a lumbar support or on top of the standard pad as a head pad. The pad can also easily be adjusted forward/back.



1400 Head Rest

#30-14-0330

Fits: 1400 Seat back

This contoured head rest is easy to mount and easy to adjust. The head rest can be adjusted up/down and forward/back for optimal positioning.

Snappy Seat Laterals

#30-10-1501M 7 inches

#30-10-1501S 5 inches

Laterals for the Snappy Seat are available in 2 sizes and is sold in pairs.



2700 Series Seat Back

#30-13-2700

The 2700 Series seat back provides added trunk support for riders with low muscle tone or balance issues. The padded seat back measures 9 inches wide by 10 inches tall and be adjusted up/down the double tube support system. Quick release cam lock levers allow the 2 removable padded lateral supports adjustments up/down, forward/back and swing in/out. Plus a strap secures lateral supports in front of rider. The 2700 Seat back can be used with a saddle seat, tractor seat, or bench seat.

Hands/Handlebars

Wrist Wraps – Now Available in 5 Sizes!


- #30-21-0160XS XSmall
- #30-21-0160S Small
- #30-21-0160M Medium
- #30-21-0160L Large
- #30-21-0160XL XLarge

Wrist wraps help maintain riders hold on hand grip. Pair of soft washable cloth wraps secure around wrist with hook and loop closure. Sizes are determined by the measurement from A: top side of wrist around hand to the B: bottom side of wrist.



Measure from point A to B for correct size.

Xtra Small	Small	Medium	Large	Xtra Large
7"	8"	10"	11"	12"



Wrist Brace Holding Mitt

- | | |
|------------------------------|---|
| #30-22-0168LXXS Left XXsmall | #30-22-0168RXXS Right XXsmall |
| #30-22-0168LXS Left Xsmall | #30-22-0168RXS Right Xsmall |
| #30-22-0168LS Left Small | #30-22-0168RS Right Small 5 ½" to 6 ½" |
| #30-22-0168LM Left Medium | #30-22-0168RM Right Medium 6 ½" to 7 1/2" |
| #30-22-0168LL Left Large | #30-22-0168RL Right Large 7 ½" to 8 1/2" |
| #30-22-0168LXL Left XLarge | #30-22-0168RXL Right XLarge 8 ½" to 9 ½" |

This is a combination of a wrist brace and holding mitt. The wrist brace holding mitt helps maintain hold on the handgrip while providing wrist stability. Durable padded material for maximum wear and comfort. Available in Left or Right and 6 sizes: XXS to XL. Measurement is wrist circumference.


Brachial Plexus Kit

- #30-23-2000 AM-9 4 inch crank
 - #30-23-2001 AM-12 and AM-16 4 inch crank
- Right, left or both sides can be altered to provide a variety of positioning. Range of motion is extended or brought to 0 degrees.



Hand Strap

- #30-20-0100
- Lightweight stretchable strap to help rider maintain hold on hand grip. Measures 30" by 2" with hook and loop closure.





Vertical handgrips

Vertical Handgrips

#30-23-0343

Fits: AM-9S, AM-9XS, AM-12 and AM-16

The vertical handgrips are designed to replace the standard handgrip on the AmTryke therapeutic tricycle for riders who are unable to pronate wrists and need help keeping hands on grips.

Dual Brake

#30-24-0100

Allows rider to control both brakes with one lever. Can be mounted either left or right side.



Dual Handbrake



Hi Rise Handlebars



MX-Freestyle Handlebars



1400 ProSeries Handlebars

Handlebars

#30- 24-0001 Hi Rise Handlebars

#30-24-0002 MX-Freestyle Handlebars

#30-24-1400 1400 ProSeries Handlebars

Pedals/Feet

AM-9 Pedal Plate

#30-30-0147

Lightweight durable plate with 1 ½" heel trap, toe and instep straps. Predrilled holes for quick and easy mounting on existing pedals.



Expanding Pedals

#30-30-0322 ½" axle Fits 1400 ProSeries and 2700 Series Models

#30-30-0623 9/16" axle Fits Snappy, AM-12 and AM-16 Models

Lightweight durable plate with 1 ½" heel trap, toe and instep straps. Plate adjusts from 7 inches to 10 inches in length to accommodate length of shoe for custom fit. The position over the pedal axis is also adjustable.



Expanding Pedals; Standard on ProSeries 1416 and 1420

Exercise Pedals

#30-30-0121 ½ " axle Fits 1400 ProSeries and 2700 Series Models

#30-30-0122 9/16" axle Fits Snappy, AM-12 and AM-16 Models

Exercise pedals provide riders that need a little assistance in keeping feet on pedals. Made with durable plastic with a texture grip and wide adjustable strap to secure riders foot in place. Replaces a standard pedal. Two axle shaft sizes available: ½ inch or 9/6 inch.



Footcups

Medium #30-31-0149M

Small #30-31-0149S


Molded from durable plastic. Medium foot cups fit children approximately 4 to 12 years of age. Small foot cups fit children approximately 1 to 3 years of age



Pedal Block

#30-31-0154

Lightweight ¾ inch wooden blocks are used to adjust for leg length differences or as an aid to help rider reach the pedal. Hook and loop securely position blocks on foot plate. Note: Using more than one pedal block will lose the heel trap on the foot plate, so foot cups would be recommended. Maximum 3 blocks per side.



AM Foot Platform

#30-32-0341

Fits: AM-12, AM-16 only

The AM foot platform provides added foot stabilization and support while functioning as a footrest. This platform allows the AM-12 and AM-16 hand / foot cycles to be used as a hand cycle for those riders that do not need the foot rotation component. The standard pedals on AM-12 and AM-16 are removed and mounted on the AM Foot Platform. The AM Foot Platform can be mounted on top or bottom of the frame and easily slide up and down the frame for foot positioning. Each foot rest can be adjusted for leg length. Note: When using the AM Foot Platform the hinged front end needs to be in the front wheel extended position.

AM-9 Pedal Leveler

#30-33-0620

Fits AM-9

The AM-9 Pedal Leveler Pulley assists the rider to maintain pedal balance. Pulley attaches to the rear of the AM-9 and cord attaches to the heel of the pedal.

Pedal Leveler Toe Pulley

#30-33-0621

Pedal Leveler Toe Pulley is an additional assist for expanding pedals on the 1400 ProSeries trykes. Pulley attaches to the front downtube of the 1400 ProSeries and cord attaches to the toe of the pedals to keep them level.

2700 Pedal Leveler

#30-33-2700

Fits: 2700 Series Models

The 2700 Pedal Leveler helps the rider maintain pedal balance. Can be used with expanding pedals only. Pulley attaches to rear of 2700 Series tryke and cord attaches to the rear of the expanding pedals.

Knee Adductor Positioning Strap

#30-40-0702S 6-9 inch circumference

#30-40-0702M 10-14 inch circumference

#30-40-0702L 15-19 inch circumference

Hook and loop strap for riders who need help with keeping their knees together so that the hamstring muscles can push down against the pedal. Available in three sizes. Sizing process starts with the circumference measurement of the thigh just above the knee.



Knee Separator

#30-41-0310 Knee Separator with 6 Inch tube

#30-41-0311 9 inch tube

Fits: AM-12 and AM-16, ProSeries 1412, 1416, and 1420

Knee separator functions as pommel or abductor. Prevents knee scissoring and keeps legs in proper alignment. Mounting bracket fits on all AM-12, AM-16 and 1400 ProSeries AmTrykes. Quick release cam lever allows easy addition of knee separator after rider is mounted. Knee separator comes standard with 6 inch mounting tube. Also available is a 9 inch mounting tube #30-41-0311.

Pull Steering Bar

#30-51-0306

Fits: Snappy, AM-12 or AM-16, and ProSeries

The pull/steering bar is for parents or therapists who wish to face the rider. It is a heavy duty one piece bar that quickly attaches to the front of the Snappy, AM-12 or AM-16.



Push Bar

Fits: Snappy, AM-12, AM-16, 1400 ProSeries 1412, 1416, 1420

#30-50-0505 2 Piece Push Bar

#30-50-0307 Heavy Duty 1 Piece Push Bar

Optimal height for use of push bar is 36 inches from handgrip to ground. The push bar is used only for assisting in movement of the AmTryke® therapeutic tricycle. It is designed to save parents or therapists from bending over when assisting a child riding the tryke. It is not intended to provide leverage for raising or lifting the tryke. #30-50-0505 is an aluminum two-sectioned bar with black plastic handgrip. #30-50-0307 is a heavy duty one piece push bar. Both are easily installed on the rear frame of the AmTryke® therapeutic tricycle.



Ordering Information

AmTryke, LLC
PO Box 5127, High Point NC 27262
1-888-AMTRYKE (268-7953)
amtryke@ambucs.org | Visit our web site: www.ambucs.org

Order Information

Please provide AmTryke with your full information when ordering to minimize any potential errors. Existing customers will be asked to verify shipping and billing information to insure proper delivery of merchandise.

Correction of Errors or Inaccuracies

The information in this catalogue may contain typographical errors or inaccuracies and may not be complete or current. We therefore reserve the right to correct any errors, inaccuracies or omissions and to change or update information on any item without prior notice (including after you have submitted your order). We apologize for any inconvenience this may cause you.

Terms

Terms are C.O.D., cashiers check, Visa, Mastercard or American Express credit cards until an open account or approval has been established. A credit application must be completed and then approved by our credit department. Returned check fee is \$35.00.

Freight Terms

Shipping and handling will be billed on your invoice. Please contact us for shipping and handling fees for each tryke model and accessory. Customers using their own shipping account numbers will be assessed a handling fee.

Returns

Merchandise may NOT be returned for credit without prior authorization from AmTryke, LLC. No returns will be accepted without prior approval. Please contact 1-888-268-7953 for return merchandise authorization number (RMA), and shipping instructions. Merchandise must be sent back within 30 days. AmTryke LLC will not be responsible for errors on size or other specifications. If you wish to return merchandise for exchange or credit, disassemble the item and pack it carefully in the original container to protect it from damage. Insure the shipment. Claims for damages or shortages must be made within five (5) days of receiving the merchandise. If you decide to cancel an order, or return an item that has already shipped or refuse an item once it arrives to you, you will be charged the freight and handling fee. A restocking fee of \$50 is charged for returned goods.

Backorders and Package Shortages

Backorders-We do not hold orders unless specifically requested. Backordered items are shipped as soon as they are received. Your order may be split into multiple invoices to allow us to fill your order in the timeliest manner possible.

Package Shortages-Claims for damages or shortages must be made within five (5) days of receiving the merchandise. Before calling to report a shortage, please make sure you have opened and thoroughly inspected ALL cartons and packing material against the packing slip that is on the outside of the box.

Shipping Methods

As a default, all orders will be shipped "best way" as determined by our shipping department to ensure prompt service and a fair price. We also offer to ship next day, second day, Saturday delivery, etc. according to your specific instructions. However, these special methods DO NOT qualify for any freight allowance.

The AmTryke Therapeutic Tricycle Program

The AmTryke Therapeutic Tricycle program is about abilities, not disabilities. The AmTryke program's goal is to provide people with disabilities with the opportunity to ride bikes. AmTryke has designed a wide range of models, sizes, and adaptations suitable for many diagnoses at low cost.

AmTryke now builds foot cycles, hand cycles, and combination hand/foot driven cycles. To insure a proper fit AmTrykes can be adapted in many ways including frame size, seat selection, trunk and leg support. This makes it ideal for children with a variety of special needs. As a therapeutic tricycle, AmTrykes help children improve motor skills, provides strength training, is fun to ride and offers the opportunity for inclusive play. For adults, AmTrykes provide a new form of mobility, exercise, and freedom to explore the community.

For therapists the use of the therapeutic tricycle will provide the clinician an age appropriate, safe, dynamic modality for improving balance, coordination, and strength. The trykes can build endurance, self confidence and provide another form of self mobility aside from the assistive device that the child may routinely use (ie. wheelchair, crutches, walker). The multiple adjustments and accessories that this unique tricycle offers makes it a valuable addition to the medical therapy

units repertoire of therapeutic devices used to further develop the functional skills of the people they serve.

For parents and families the use of the AmTryke fulfills the need of every child to have a bike. All children deserve a chance to have a bike just like their siblings and friends. Given a safe, stable platform, many riders not previously considered can ride just like their peers. AmTrykes are designed to look like bikes, not medical equipment, and to be age appropriate.

AMBUCS has provided over 12,000 AmTrykes to riders across America. A grass roots charity, AMBUCS chapters fundraise and provide the equipment for those with financial need. Riders must be evaluated by a physical or occupational therapist to be eligible.

The AmTryke Road Trip brings the equipment to rehabs, clinics, and therapy centers. Trainings or inservices are offered to teach therapists about the AmTryke program. These are often combined with bike days so families can bring their riders and receive free evaluations for AmTrykes.

Contact the AMBUCS Resource Center to schedule a visit from The AmTryke Road Show today!